

COMPREHENSIVE VALIDATION PACKAGE

ATL Applications
INVENTORY SHEET

WORK ORDER # 0908456C

	Page Nos.	
	From	To
1. Work Order Cover Page & Laboratory Narrative & Table	1	4
2. Sample Results and Raw Data (Organized By Sample)	5	8
a. ATL Sample Results Form		
b. Target Compound Raw Data		
-Internal Standard Area and Retention Time Summary (If Applicable)		
-Surrogate Recovery Summary (If Applicable)		
-Chromatogram(s) and Ion Profiles (If Applicable)		
3. QC Results and Raw Data		
a. Method Blank (Results + Raw Data)	-	-
b. Surrogate Recovery Summary Form (If Applicable)	-	-
c. Internal Standard Summary Form (If Applicable)	-	-
d. Duplicate Results Summary Sheet	-	-
e. Matrix Spike/Matrix Spike Duplicate (Results + Raw Data)	-	-
f. Initial Calibration Data (Summary Sheet + Raw Data)		
g. MDL Study (If Applicable)	-	-
h. Continuing Calibration Verification Data		
i. Second Source LCS (Summary + Raw Data)	-	-
j. Extraction Logs	-	-
k. Instrument Run Logs/Software Verification	9	14
l. GC/MS Tune (Results + Raw Data)	-	-
4. Shipping/Receiving Documents:		
a. Login Receipt Summary Sheet	15	16
b. Chain-of-Custody Records	17	18
c. Sample Log-In Sheet	19	20
d. Misc. Shipping/Receiving Records (list individual records)		
<u>Sample Receipt Discrepancy Report</u>	21	23
5. Other Records (describe or list)		
a. <u>Manual Spectral Defense</u>	-	-
b. <u>Manual Intergrations</u>	-	-
c. <u>Manual Calculations</u>	-	-
d. <u>Canister Dilution Factors</u>	-	-
e. <u>Laboratory Corrective Action Request</u>	-	-
f. <u>CAS Number Reference</u>	24	25
g. <u>Variance Table</u>	-	-
h. <u>Canister Certification</u>	-	-
i. <u>Data Review Check Sheet</u>	26	26

Completed by:

Kara McKiernan

(Signature)

Kara McKiernan/ Document Control

(Print Name & Title)

09/17/09

(Date)

WORK ORDER #: 0908456C

Work Order Summary

CLIENT:	Mr. Taeko Minegishi Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494	BILL TO:	Accounts Payable Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494
PHONE:	800-825-5343	P.O. #	16512
FAX:	781-247-4305	PROJECT #	16512
DATE RECEIVED:	08/21/2009	CONTACT:	Ausha Scott
DATE COMPLETED:	09/16/2009		

<u>FRACTION #</u>	<u>NAME</u>	<u>TEST</u>
41A	100600	ATL Applications
41AA	100600 Lab Duplicate	ATL Applications
42A	100601	ATL Applications
43A	100602	ATL Applications
44A	100603	ATL Applications
45A	100604	ATL Applications
46A	100605	ATL Applications
47A	100171	ATL Applications
48A	100172	ATL Applications
49A	100173	ATL Applications
50A	100174	ATL Applications
50AA	100174 Lab Duplicate	ATL Applications
51A	100175	ATL Applications
52A	100176	ATL Applications
53A	100231	ATL Applications
54A	100232	ATL Applications
55A	100233	ATL Applications

Continued on next page


WORK ORDER #: 0908456C

Work Order Summary

CLIENT:	Mr. Taeko Minegishi Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494	BILL TO:	Accounts Payable Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494
PHONE:	800-825-5343	P.O. #	16512
FAX:	781-247-4305	PROJECT #	16512
DATE RECEIVED:	08/21/2009	CONTACT:	Ausha Scott
DATE COMPLETED:	09/16/2009		

<u>FRACTION #</u>	<u>NAME</u>	<u>TEST</u>
56A	100234	ATL Applications
57A	100235	ATL Applications
58A	100236	ATL Applications
59A	Method Blank	ATL Applications
59B	Method Blank	ATL Applications
60A	CCV	ATL Applications

CERTIFIED BY:



Laboratory Director

DATE: 09/16/09

This report shall not be reproduced, except in full, without the written approval of Air Toxics Ltd.

180 BLUE RAVINE ROAD, SUITE B FOLSOM, CA - 95630
(916) 985-1000 . (800) 985-5955 . FAX (916) 985-1020

LABORATORY NARRATIVE
Hydrogen Sulfide by Radiello 170
Environmental Health & Engineering, Inc.
Workorder# 0908456C

Eighteen Radiello 170 (H₂S) samples were received on August 21, 2009. The procedure involves adsorption of H₂S by zinc acetate to form zinc sulfide. The sulfide is then recovered by extraction with water and addition of ferric chloride in a strongly acidic solution to produce methylene blue. Methylene blue absorbance is then measured at 665 nm using a spectrophotometer. Results are reported in uG and uG/m³.

Sampling rate of 69 mL/min for H₂S was provided by the manufacturer.

Receiving Notes

A Temperature Blank was not included with the shipment. Temperature was measured on a representative sample and was not within 4±2 °C. Coolant in the form of blue ice was present. Analysis proceeded.

Sample collection dates were not provided on the Chain of Custody for all samples. The client was contacted and a dates were provided.

Analytical Notes

Results were calculated based on 25 deg C without temperature correction. The actual exposure time was used to calculate sample concentrations and reporting limits.

An exposure time of 20000 minutes was used for the QC samples.

All media used for the sampling were supplied by the client. Blank subtraction was not performed on the sample results since the media used for Method Blanks may be from a different lot than the media used for the samples.

Definition of Data Qualifying Flags

Eight qualifiers may have been used on the data analysis sheets and indicate as follows:

B - Compound present in laboratory blank greater than reporting limit.

J - Estimated value.

E - Exceeds instrument calibration range.

S - Saturated peak.

Q - Exceeds quality control limits.

U - Compound analyzed for but not detected above the detection limit.

M - Reported value may be biased due to apparent matrix interferences.

N - The identification is based on presumptive evidence.

File extensions may have been used on the data analysis sheets and indicates as follows:

a-File was requantified

b-File was quantified by a second column and detector

r1-File was requantified for the purpose of reissue

Sample Results and Raw Data

AIR TOXICS LTD.

ATL Application # 59 for RAD 170 (Hydrogen Sulfide)

Spectrophotometer

Field Sample ID.	Lab Sample ID.	Collection Date	Analysis Date	Dilution Factor	Reporting Limit (ug)	Reporting Limit (ug/m3)	Amount (ug)	Amount (ug/m3)
100600	0908456C-41A	8/19/2009	8/25/2009	1.00	0.80	0.55	5.5	3.8
100600 Duplicate	0908456C-41AA	8/19/2009	8/25/2009	1.00	0.80	0.55	5.5	3.8
100601	0908456C-42A	8/19/2009	8/25/2009	1.00	0.80	0.55	1.4	0.94
100602	0908456C-43A	8/19/2009	8/25/2009	1.00	0.80	0.55	3.9	2.7
100603	0908456C-44A	8/19/2009	8/25/2009	1.00	0.80	0.55	3.8	2.6
100604	0908456C-45A	8/19/2009	8/25/2009	1.00	0.80	0.55	ND	ND
100605	0908456C-46A	8/19/2009	8/25/2009	1.00	0.80	0.55	ND	ND
100171	0908456C-47A	8/18/2009	8/25/2009	1.00	0.80	0.55	1.4	0.96
100172	0908456C-48A	8/18/2009	8/25/2009	1.00	0.80	0.55	1.6	1.1
100173	0908456C-49A	8/18/2009	8/25/2009	1.00	0.80	0.55	2.8	1.9
100174	0908456C-50A	8/18/2009	8/25/2009	1.00	0.80	0.55	1.8	1.2
100174 Duplicate	0908456C-50AA	8/18/2009	8/25/2009	1.00	0.80	0.55	1.8	1.2
100175	0908456C-51A	8/18/2009	8/25/2009	1.00	0.80	0.55	1.5	1.0
100176	0908456C-52A	8/18/2009	8/25/2009	1.00	0.80	0.55	ND	ND
100231	0908456C-53A	8/20/2009	8/25/2009	1.00	0.80	0.55	1.9	1.3
100232	0908456C-54A	8/20/2009	8/25/2009	1.00	0.80	0.55	1.8	1.3
100233	0908456C-55A	8/20/2009	8/25/2009	1.00	0.80	0.55	1.6	1.1
100234	0908456C-56A	8/20/2009	8/25/2009	1.00	0.80	0.55	ND	ND
100235	0908456C-57A	8/20/2009	8/25/2009	1.00	0.80	0.55	ND	ND
100236	0908456C-58A	8/20/2009	8/25/2009	1.00	0.80	0.55	ND	ND
Method Blank	0908456C-59A	NA	8/25/2009	1.00	0.80	0.55	ND	ND
Method Blank	0908456C-59B	NA	8/25/2009	1.00	0.80	0.55	ND	ND
CCV	0908456C-60A	NA	8/25/2009	1.00	0.80	0.55	%Rec 100	

COMMENTS: 1. NA=Not Applicable
 2. ND=Not Detected
 3. Exposure time of 20000 minutes was assumed for the QC samples and samples 100604, 100605, 100176 and 100236.
 4. Background subtraction not performed.

Hydrogen Sulfide Radiello Calculation Worksheet

Workorder #: 0908456C

Sampling Rate (mg/ppb-min)

0.096 Typically 0.096 for H2S

Sampling T (deg C)

25 Typically 25

Volume (ml)

10.5 Typically 10.5 for H2S

Date of Analysis:

8/25/2009

Corrected Q

0.096 Takes into account temp

(Abs-Y-int)xDF

Conc(ug/ml)xVol (ml)

conc (ug sulfide) *MW H2S

Q includes conversion from Sulfide to H2S

Conc (ug) x 1000

ppbx mw

T Corrected, no Blank correction

LabSampleID	Client	Date of Collection	Abs	Duration (min)	DF	Conc (ug/ml) of sulfide	Conc (ug) of sulfide	Conc (ug) of H2S	Conc (ppb) of H2S	Conc (ug/m3) of H2S
41A	100800	8/19/2009	0.539	19850	1.00	0.493674332	5.18358049	5.508786775	2.720	3.792
41AA	100800 Duplicate	8/19/2009	0.538	19850	1.00	0.492717189	5.173530481	5.498106201	2.715	3.784
42A	100801	8/19/2009	0.151	19850	1.00	0.122302558	1.284176859	1.36474324	0.674	0.939
43A	100802	8/19/2009	0.384	19850	1.00	0.345317051	3.62582904	3.853305436	1.903	2.652
44A	100803	8/19/2009	0.376	19850	1.00	0.337659901	3.545428965	3.76786124	1.861	2.593
45A	100804	8/19/2009	0.015	19963	1.00	-0.007868992	-0.082624414	-0.087808085	-0.043	-0.060
46A	100805	8/19/2009	0.013	19963	1.00	-0.009783279	-0.102724433	-0.109169134	-0.054	-0.075
47A	100171	8/18/2009	0.154	19963	1.00	0.125173989	1.314326887	1.396784813	0.686	0.956
48A	100172	8/18/2009	0.173	19963	1.00	0.14335972	1.505277065	1.599714778	0.785	1.085
49A	100173	8/18/2009	0.287	19963	1.00	0.252474108	2.650978132	2.817294565	1.383	1.928
50A	100174	8/18/2009	0.188	19963	1.00	0.157716877	1.656027205	1.759922644	0.864	1.205
50AA	100174 Duplicate	8/18/2009	0.188	19963	1.00	0.157716877	1.656027205	1.759922644	0.864	1.205
51A	100175	8/18/2009	0.163	19963	1.00	0.133788283	1.404776971	1.492905533	0.733	1.022
52A	100176	8/18/2009	0.014	19963	1.00	-0.008826136	-0.092674424	-0.098488609	-0.048	-0.067
53A	100231	8/20/2009	0.198	19712	1.00	0.167288314	1.756527299	1.866727889	0.928	1.294
54A	100232	8/20/2009	0.193	19712	1.00	0.162502595	1.70627752	1.813325267	0.902	1.257
55A	100233	8/20/2009	0.171	19712	1.00	0.141445433	1.485177046	1.578353729	0.785	1.094
56A	100234	8/20/2009	0.033	19712	1.00	0.009359596	0.098275754	0.104441355	0.052	0.072
57A	100235	8/20/2009	0.095	19712	1.00	0.068702508	0.721376335	0.766633871	0.381	0.531
58A	100236	8/20/2009	0.014	19963	1.00	-0.008826136	-0.092674424	-0.098488609	-0.048	-0.067
59A	Method Blank	NA	0.018	19963	1.00	-0.0222226148	-0.233374555	-0.248015952	#DIV/0!	#DIV/0!
59B	Method Blank	NA	0.016	19963	1.00	-0.004997561	-0.052474386	-0.055766512	-0.027	-0.038
60A	CCV	NA	0.618	19963	1.00	0.569288689	5.97753123	6.352548156	-0.038	-0.053
				QC Duration		CCV Spike Amt				
				19963		0.572			3.119	4.348

Low Point:DF

RL(ug/ml)XVol (ml)

RL (ug sulfide) * MW H2S
MW SulfideQ includes conversion from
Sulfide to H2SRL (ug) x 1000
Q x DurationppbX mw
24.45

Calibration Data

Calibration Date
8/25/2009 Linear Regression

RL(ug/ml) of sulfide	RL (ug) of sulfide	Conc (ug) of H2S	RL (ppb) of H2S	RL (ug/m3)	T Corrected, no Blank correction			ug/ml of sulfide	absorbance	Slope Y-int R2
					Result (ug) H2S	Result (ug/m3) H2S	%Rec			
0.072	0.752	0.798966249	0.39	0.55	5.508786725	3.791790299		0	0	1.044775146
0.072	0.752	0.798966249	0.39	0.55	5.498106201	3.784438715		0.0716	0.088	0.023221327
0.072	0.752	0.798966249	0.39	0.55	1.36474324	0.939375662		0.143	0.169	0.999419909
0.072	0.752	0.798966249	0.39	0.55	3.853305436	2.652294761		0.286	0.328	
0.072	0.752	0.798966249	0.39	0.55	3.76786124	2.593482088		0.572	0.637	
0.072	0.752	0.798966249	0.39	0.55	ND	ND		1.145	1.211	
0.072	0.752	0.798966249	0.39	0.55	1.396784813	0.955988265				
0.072	0.752	0.798966249	0.39	0.55	1.599714778	1.094877708				
0.072	0.752	0.798966249	0.39	0.55	2.817294565	1.928214365				
0.072	0.752	0.798966249	0.39	0.55	1.759922644	1.204527268				
0.072	0.752	0.798966249	0.39	0.55	1.759922644	1.204527268				
0.072	0.752	0.798966249	0.39	0.55	1.492909533	1.021778001				
0.072	0.752	0.798966249	0.40	0.55	1.866727889	1.293895459				
0.072	0.752	0.798966249	0.40	0.55	1.813325267	1.256880203				
0.072	0.752	0.798966249	0.40	0.55	1.578353729	1.094013078				
0.072	0.752	0.798966249	0.40	0.55	ND	ND				
0.072	0.752	0.798966249	0.39	0.55	ND	ND				
0.072	0.752	0.798966249	#DNV/OI	#DNV/OI	ND	#DNV/OI				
0.072	0.752	0.798966249	#DNV/OI	#DNV/OI	ND	#DNV/OI				
0.072	0.752	0.798966249	0.39	0.55	ND	ND				
0.072	0.752	0.798966249	0.39	0.55	6.352548156	4.347814661	%Rec			
0.072	0.752	0.798966249	0.39	0.55			100			

QC Results and Raw Data

Spectrophotometer Logbook

@Air Toxics Ltd.

Log Book #: 1564Work Order: 0908456 UDDate: 8/25/09Method: Rad 170Analyst: A. ToyamaWavelength: 665

Prep. Notes:

Standard ID	Concentration	ABS
1858-20-E	0.0716 ug/mL	0.088
-D	0.143	0.169
-C	0.286	0.328
-B	0.572	0.637
-A	1.145	1.211

$$r = 0.9994199$$

$$m = 1.044775$$

$$b = 0.023221$$

Fraction	Dilution	ABS	Sample ID	Sample Volume
41A	1.00	0.539	100600	10.5 mL
42A		0.151	601	
43A		0.384	602	
44A		0.376	603	
45A		0.015	604	
46A		0.013	605	
47A		0.154	171	
48A		0.173	172	
49A		0.287	173	
50A		0.188	174	
51A		0.163	175	
52A		0.014	176	
53A		0.198	231	
54A		0.193	232	
55A		0.171	233	

Notes: Code 170 Lot 09075 Exp 04/010 used for Blanks. Sample lots unknown.

Log Book #: 1564

Prep. Notes: cont. from page 28

Standard ID	Concentration	ABS
1858-20-E	0.0716 $\mu\text{g/mL}$	0.088
-D	0.143	0.169
-C	0.286	0.328
-B	0.572	0.637
-A	1.145	1.211

r = _____
m = _____
b = _____

[illegible]

Notes: LCS/CCV prepared at 0.572 g/ml

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1858

Standard ID: 1858-18

Project: Ferric chloride ~~amine~~ 8/24/09 ACT RAD70

Analyst: A. Toyama

Solvent: DI H₂O

Solvent Lot #: NA

Preparation Date: 8/25/09

Expiration Date: 8/25/09

Procedure/Comments: Dissolve 25 g of ferric chloride hexahydrate (Located ER2C
Lot 73297 MS) in ¹⁰⁰~~125~~ mL of DI H₂O.
8/24/09 ACT

8/24/09
ACT

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1858

Standard ID: 1858-19

Project: Ferric chloride - amine

Analyst: A. Toxams

Preparation Date: 8/25/09

Expiration Date: 8/25/09

Solvent: DI H₂O

Solvent Lot #: NA

Procedure/Comments: Mix 10 mL of ferric chloride solution with 50 mL of amine solution.

8/24/09
ACT

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1858

Standard ID: 1858-20

Project: Calibration Solution RAD 170

Analyst: A. Toyama

Preparation Date: 8/25/09

Expiration Date: 8/25/09

Solvent: DI H₂O

Solvent Lot #: NA

Procedure/Comments: _____

Solution A: 2 mL of Code RAD 171 (1476-984, Exp 8/6/10, Located ER13) with 98 mL
DI H₂O = 1.145 $\mu\text{g/mL}$ sulfide ions

Solution B: 25 mL of Solution A with 25 mL of DI H₂O = 0.572 $\mu\text{g/mL}$

Solution C: 1.25 mL of Solution A with 3.75 mL of DI H₂O = 0.286 $\mu\text{g/mL}$

Solution D: 0.625 mL of Solution A with 4.375 mL of DI H₂O = 0.143 $\mu\text{g/mL}$

Solution E: 0.375 ^{mL} of Solution A with 5.625 mL of DI H₂O = 0.0716 $\mu\text{g/mL}$
- (remove 1.0 mL of final solution) 8/25/09
AK

8/24/09
AK

Shipping/ Receiving Documents

**180 Blue Ravine Road, Suite B
Folsom, CA 95630**

**Phone (916) 985-1000 FAX (916) 985-1020
Hours 8:00 A.M. to 6:00 P.M. Pacific**

COMPANY: Environmental Health & Engineering, Inc.
ATTENTION: Mr. Taeko Minegishi
FAX #: 781-247-4305
FROM: Sample Receiving
Workorder #: 0908456C
of pages (Including Cover): 4

9/17/2009

Thank you for selecting Air Toxics Ltd. We have received your samples and have found no discrepancies.
In order to expedite analysis and reporting, please review the attached information for accuracy.

Corrections can be faxed to **Ausha Scott at 916-985-1020.**

ATL will proceed with the analysis as specified on the Chain of Custody and Sample Login page.

FROM: Environmental Health and Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494-2725

TO: AIR TOXICS

Please send invoices to ATTN: Accounts Payable
Please send reports to ATTN: Data Coordinator

In all correspondence regarding this matter, please refer to EH&E Project # 16512

The cost of this analysis will be covered by EH&E Purchase Order # 16512

For EH & E Data Coordinator - URGENT DATA ☒

	SAMPLE ID	SAMPLE TYPE	ANALYTICAL METHOD/NUMBER	OTHER: Time/Date/Vol.
33A	100348	AIR PASSIVE	H ₂ S ANALYSIS	φ
34A	100698			12D 18H 52M
35A	100699			
36A	100700			
37A	100701			
38A	100702			
39A	100703			φ
40A	100599			13D 18H 50M
41A	100600			
42A	100601			
43A	100602			
44A	100603			
45A	100604			φ
46A	100605			φ
47A	100171			13D 20H 43M
48A	100172			

Special Instructions:

☒ Standard turn around time


☐ Rush by _____ date/time

☐ Fax results 781-247-4305

☐ RETURN SAMPLES

☒ Electronic transfer - datacoordinator@eh&e.com

☒ Additional report recipient mragala@eh&e.com



Each signatory please return one copy of this form to the above address

Relinquished by: [Signature] of Environmental Health & Engineering, Inc. Date: 8/21/09

Received by: [Signature] 0850 of (company name) ATL Date: 8/21/09

Relinquished by: _____ of (company name) _____ Date: _____

Received by: _____ of (company name) _____ Date: _____

Relinquished by: _____ of (company name) _____ Date: _____

Received by: _____ of (company name) _____ Date: _____

Lab Data

Received by: _____ of Environmental Health & Engineering, Inc. Date: _____

Page 3 of 4

FROM: Environmental Health and Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494-2725

TO: AIR TOXICS

Please send invoices to ATTN: Accounts Payable
Please send reports to ATTN: Data Coordinator

In all correspondence regarding this matter, please refer to EH&E Project # 16512

The cost of this analysis will be covered by EH&E Purchase Order # 16512

For EH & E Data Coordinator - URGENT DATA ☒

SAMPLE ID	SAMPLE TYPE	ANALYTICAL METHOD/NUMBER	OTHER: Time/Date/Vol.
49A 100173	AIR/PASSIVE	H ₂ S ANALYSIS	130 20H 43M
50A 100174			
51A 100175			
52A 100176			
53A 100231			130 16H 32M
54A 100232			
55A 100233			
56A 100234			
57A 100235			
58A 100236			

Special Instructions:

- ☒ Standard turn around time ☐ Rush by _____ date/time ☐ Other _____
☐ Fax results 781-247-4305 ☒ Electronic transfer - datacoordinator@eh&e.com
☐ RETURN SAMPLES ☒ Additional report recipient wfrayala@eh&e.com

Each signatory please return one copy of this form to the above address

Relinquished by: [Signature] of Environmental Health & Engineering, Inc. Date: 8/21/09

Received by: AR 0830 of (company name) ATC Date: 8/21/09

Relinquished by: _____ of (company name) _____ Date: _____

Received by: _____ of (company name) _____ Date: _____

Relinquished by: _____ of (company name) _____ Date: _____

Received by: _____ of (company name) _____ Date: _____

Lab Data

Received by: _____ of Environmental Health & Engineering, Inc. Date: _____

Page 4 of 4

SAMPLE RECEIPT SUMMARY

WORKORDER 0908456C

Client

Mr. Taeko Minegishi
Environmental Health &
Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494

Phone

800-825-5343

Fax

781-247-4305

Date Promised: 09/01/09 11:59 pm

Date Completed: 9/16/09

Date Received: 8/21/09

PO#: 16512

Project#: 16512

Sales Rep: TL

Total \$: \$ 1,000.00

Logged By: MW

<u>Fraction</u>	<u>Sample #</u>	<u>Analysis</u>	<u>Collected</u>	<u>Amount\$</u>
41A	100600	ATL Applications	8/19/2009	\$50.00
41AA	100600 Lab Duplicate	ATL Applications	8/19/2009	\$0.00
42A	100601	ATL Applications	8/19/2009	\$50.00
43A	100602	ATL Applications	8/19/2009	\$50.00
44A	100603	ATL Applications	8/19/2009	\$50.00
45A	100604	ATL Applications	8/19/2009	\$50.00
46A	100605	ATL Applications	8/19/2009	\$50.00
47A	100171	ATL Applications	8/18/2009	\$50.00
48A	100172	ATL Applications	8/18/2009	\$50.00
49A	100173	ATL Applications	8/18/2009	\$50.00
50A	100174	ATL Applications	8/18/2009	\$50.00
50AA	100174 Lab Duplicate	ATL Applications	8/18/2009	\$0.00
51A	100175	ATL Applications	8/18/2009	\$50.00
52A	100176	ATL Applications	8/18/2009	\$50.00
53A	100231	ATL Applications	8/20/2009	\$50.00
54A	100232	ATL Applications	8/20/2009	\$50.00
55A	100233	ATL Applications	8/20/2009	\$50.00
56A	100234	ATL Applications	8/20/2009	\$50.00
57A	100235	ATL Applications	8/20/2009	\$50.00
58A	100236	ATL Applications	8/20/2009	\$50.00

Note: Samples received after 3 P.M. PST are considered to be received on the following work day.
Atlas Project Name/Profile#: CPSC Indoor Air Monitoring/13297

BILL TO: Accounts Payable
Environmental Health & Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494

Analysis Code: Other GC

TERMS:

Reporting Method: ATL Application #59 H2S-Radiello 170

180 BLUE RAVINE ROAD, SUITE B FOLSOM, CA - 95630
(916) 985-1000 . (800) 985-5955 . FAX (916) 985-1020

SAMPLE RECEIPT SUMMARY Continued

Client	Phone	Date Promised: 09/01/09 11:59 pm
Mr. Taeko Minegishi	800-825-5343	Date Completed: 9/16/09
Environmental Health &		Date Received: 8/21/09
Engineering, Inc.	Fax	PO#: 16512
117 Fourth Avenue	781-247-4305	Project#: 16512
Needham, MA 02494		
Sales Rep: TL		Total \$: \$ 1,000.00
		Logged By: MW

<u>Fraction</u>	<u>Sample #</u>	<u>Analysis</u>	<u>Collected</u>	<u>Amount\$</u>
59A	Method Blank	ATL Applications	NA	\$0.00
59B	Method Blank	ATL Applications	NA	\$0.00
60A	CCV	ATL Applications	NA	\$0.00
Misc. Charges eCVP (20) @ \$5.00 each.				\$100.00

Note: Samples received after 3 P.M. PST are considered to be received on the following work day.
Atlas Project Name/Profile#: CPSC Indoor Air Monitoring/13297

BILL TO: Accounts Payable
Environmental Health & Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494

Analysis Code: Other GC

TERMS:

Reporting Method: ATL Application #59 H2S-Radiello 170

180 BLUE RAVINE ROAD, SUITE B FOLSOM, CA - 95630
(916) 985-1000 . (800) 985-5955 . FAX (916) 985-1020

Sample Discrepancy Report

Identification

Initiated By: MW Project ID: 13297 PM: BL Date: 8/21/2009 Discrepancy Type: ☐ 1. ☒ 2. ☐ 3.

Workorder(s) affected: 0908456 Sample(s) affected: all

1. Sample Receipt Discrepancies

Narration Not Required:

- 1.1. ☐ Sample container (cartridge/tube/VOA vial) was received broken, however sample was intact.
- 1.2. ☐ No brass cap on canister.
- 1.3. ☐ Date of Collection noted on first sample, but no arrow down to indicate all samples.

Notify Lab for further determination:

- 1.4. ☐ Tedlar bag received with minimal volume.

Initials: _____ Date: _____

Narration Required in Lab Narrative and Sample Confirmation:

- 1.5. ☐ COC was not filled out in ink.
- 1.6. ☐ COC improperly relinquished / received.
- 1.7. ☐ Sample tags / can numbers do not match the COC.
- 1.8. ☐ Sample date ☐ error / ☐ missing on COC but noted on sample tag (check one).
- 1.9. ☐ Custody Seal on the outside of the container was ☐ broken / ☐ improperly placed (check one).
- 1.10. ☐ ID-none on the sample Tag/Blank
- 1.11. ☐ Other (describe below).

Describe the Discrepancy:

2. Sample Receipt/Screening Discrepancies requiring PM notification

Document on Cover Page of Sample Receipt Confirmation and in Receiving Notes of Lab Narrative

If Section II. is filled out PM must be notified within 24 hrs of initiation

- 2.1. ☐ COC was not received with samples.
- 2.2. ☐ Analysis method(s) is ☐ not specified / ☐ incorrectly specified (check one) on the COC.
- 2.3. ☐ Incorrect sampling media / container for analysis requested.
- 2.4. ☐ Number of samples on the COC does not match the number of samples that were received.
- 2.5. ☐ Samples were received expired.
- 2.6. ☒ Sampling date (time for sulfur) is not documented for ☐ some / ☒ any samples (check one).
- 2.7. ☐ Sample received with amount of H₂O in the Tedlar Bag.
- 2.8. ☐ Sample cannot be analyzed. Container was ☐ received broken / ☐ leaking / ☐ flat / ☐ defective.
- 2.9. ☐ Tedlar bag / canister received emitting a strong odor; Sample ☐ can / ☐ cannot (check one) be analyzed.
- 2.10. ☐ Tedlar Bag for Sulfur analysis has metal fitting.
- 2.11. ☐ Environmental Supply Company valves
- 2.12. ☐ Sorbent samples-sampling volume was not provided
- 2.13. ☐ Flow controller used – canister samples received at ambient or under pressure.
- 2.14. ☐ Canister was at ambient pressure at time of pressurization and (check all that apply):
 - ☐ Canister failed leak check on two manifolds,
 - ☐ Canister valve was open,
 - ☐ Brass nut was loose/not present.
 - ☐ Sample can be analyzed
 - ☐ Cannot be analyzed
- 2.15. ☐ Canister sample received with a vacuum difference >5.0"Hg between the receipt vac. And the final vac. reported on the COC, indicating loss of vacuum.
- 2.16. ☐ Canister sample received at >15"Hg (not identified as a Trip/Field Blank).
- 2.17. ☐ Canister Trip Blank received at low vacuum (< 25"Hg).
- 2.18. ☒ Sorbent Sample received outside method required temperature of 2°C to 6°C; ☐ ice / ☒ blue ice (check one) was present. A temp. Blank ☐ was / ☒ was not present (check one).
- 2.19. ☐ Other (describe below)

Initials: _____

Date: _____

Notify Receiving: ☐

Notify PM: ☐

Describe the Discrepancy: samples rec'd at 8C

3. Lab Discrepancies requiring Team Leader/PM notification

Document in Analytical Notes of Lab Narrative

If Section III. is filled out PM must be notified within 24 hrs of initiation

- | | |
|--|--|
| 3.1. <input type="checkbox"/> Tedlar Bag found to be leaking at the time of analysis; sample <input type="checkbox"/> can / <input type="checkbox"/> cannot (check one) be analyzed. | 3.6. <input type="checkbox"/> Sample loss due to instrument malfunction / broken glassware. |
| 3.2. <input type="checkbox"/> Tedlar Bag found to be flat/low volume; sample cannot be analyzed. | 3.7. <input type="checkbox"/> Low/high surrogate recoveries noted in QC/sample(s) for extractable samples. |
| 3.3. <input type="checkbox"/> Sulfur samples received with insufficient time to analyze prior to expiration. | 3.8. <input type="checkbox"/> Reporting Limit was raised. |
| 3.4. <input type="checkbox"/> Canister found to be leaking at the time of analysis. | 3.9. <input type="checkbox"/> Post weight > Pre weight in field/lab Blank for PM10/TSP samples. |
| 3.5. <input type="checkbox"/> VOST tube saturated; bag dilution necessary. | 3.10. <input type="checkbox"/> Other (describe below). |

Initials: _____ Date: _____ Notify Receiving: ☐ Notify PM: ☐

Team Lead Initials: _____ Date: _____

Describe the Discrepancy: _____

How Does this Affect Client: _____

Project Manager Use Only

Project Manager Notification Complete

☒ Section 2 Complete

☐ Section 3

Action:

- ☐ It is not necessary to notify the client. Narrate the discrepancy in Receiving Notes/Analytical Notes of Lab Narrative.

PM Initials: _____ Date: _____

- ☒ Client notification required. See attached client contact / email, or comments below:

Client Notification:

PM Initials: BL Person notified: David Shore

Date: 8/21/2009

- ☐ Waiting for Client Reply

Comments: **Proceed and narrate temperature discrepancy. See table for time of collection.**

☐ Notify Lab Name: _____ Date: _____ Notify Receiving: ☒

- ☐ Additional notifications attached.

Additional Comments:

Other Records

Method : ATL Application #59 H2S-Radiello 170

CAS Number	Compound	Rpt. Limit (ug)
7783-06-4	Hydrogen Sulfide	1.2

DATA REVIEW CHECKLIST

Work Order #:

0908456C

A₁ A₂ R T M Q

- ☐ ☐ ☒ ☐ ☒ ☐ Analysis/Reporting vs. Project Profile/SOP requirements checked (i.e. 100% Dups, J-Flag to MDL, etc)
☐ ☐ ☒ ☐ ☒ ☐ The final report has the correct reporting list, special units, and header info.
☐ ☐ ☒ ☐ ☒ ☐ Lab Narrative is correct (proper method & description/Receiving & Analytical notes correct)
☐ ☐ ☒ ☐ ☒ ☐ Sample Discrepancy Report (SDR) is completed

- ☐ ☐ ☐ ☒ ☒ ☐ Corrective Action issued - # _____
☐ ☐ ☐ ☐ ☒ ☐ Unusual circumstances have been documented in the notes section below

LUMEN validation report present and initialed

CIRCLE (YES / NO)

- ☐ ☐ ☒ ☐ ☒ ☐ Lab Blank, CCV, LCS and DUP met QC criteria
☐ ☐ ☒ ☐ ☒ ☐ Hold time is met for all samples
☐ ☐ ☒ ☐ ☒ ☐ Appropriate data qualifier flags are applied
☐ ☐ ☒ ☐ ☒ ☐ Manual integrations for samples and QC are properly documented
☐ ☐ ☒ ☐ ☒ ☐ Samples analyzed within the project or method specific clock
☐ ☐ ☒ ☐ ☒ ☐ Retention times have been verified
☐ ☐ ☒ ☐ ☒ ☐ Appropriate ICAL(s) included
☐ ☐ ☒ ☐ ☒ ☐ At least one result per sample is verified against the target quant sheets/raw data
- ☐ ☐ ☒ ☐ ☐ ☐ Dilution factor correctly calculated (sample load volume, syringe and bag dilutions, can pressurization(s))
☐ ☐ ☒ ☐ ☐ ☐ Correct amount of sample analyzed (i.e. sample not over-diluted)
☐ ☐ ☒ ☐ ☐ ☐ Spectra verified - documentation of spectral defense included (Section 5A of eCVP pkg)
☐ ☐ ☒ ☐ ☐ ☐ TICs resemble reference spectra
☐ ☐ ☒ ☐ ☐ ☐ TICs between duplicate samples are consistent
☐ ☐ ☒ ☐ ☒ ☐ Checked samples for trends (i.e. Influent vs. Effluent, Field Dups, Field/Trip Blank, etc.)
☐ ☐ ☒ ☐ ☐ ☐ Data for multiple analyses of sample(s) has been evaluated for comparability of results
☐ ☐ ☒ ☐ ☒ ☐ Special units for all samples in the final report are correctly calculated
☐ ☐ ☒ ☐ ☒ ☐ Manually entered results checked (i.e. TPH/NMOC)
☐ ☐ ☒ ☐ ☐ ☐ Chain of Custody verified for any special comments (i.e. different compounds/RLs, action levels)
☐ ☐ ☒ ☐ ☐ ☐ Chain of Custody scanned correctly
☐ ☐ ☒ ☐ ☐ ☐ Verify sample id's vs. chain of custody
☐ ☐ ☒ ☐ ☐ ☐ Date MDL(s) performed per instrument(s) _____
☐ ☐ ☒ ☐ ☐ ☐ Samples pressurized w/ appropriate gas (N₂ or He) ☒ Other (i.e. Tedlar bag, cartridge, sorbent)
☐ ☐ ☒ ☐ ☐ ☐ Final pressure consistent with canister size (6L vs. 1L)
☐ ☐ ☒ ☐ ☐ ☐ Verify receipt pressures
☐ ☐ ☒ ☐ ☐ ☐ Verify canister ID #'s
☐ ☐ ☒ ☐ ☒ ☐ Final invoice amount correct (adjusted for TAT, Penalties, Re-issue Charges etc.)
☐ ☐ ☒ ☐ ☒ ☐ MDL date(s) present for all instruments utilized
☐ ☐ ☒ ☐ ☒ ☐ Client LUMEN report reviewed for accuracy and completeness

Notes: (to include: noting samples with QA/QC problems, Blanks with positive hits, narratives, etc.)

A/R:

Dup. 41A, 50A

M/Q:

A₁/A₂

(Analytical Review/Date)

A₁:

R/T

(Reporting Review/Date)

R:

M

(Management Review/Date)

M:

Q

(QA Review/Date)

A₂:

T:

Note (1): Please check all the appropriate boxes. Indicate "NA" for any statement that does not apply.

Rev. 02/20/09

Note (2): Management reviewer and reporting reviewer must be separate individuals.