

COMPREHENSIVE VALIDATION PACKAGE


ATL Applications

INVENTORY SHEET

WORK ORDER # 0908455B

| | Page Nos. | |
|--|-----------|----|
| | From | To |
| 1. Work Order Cover Page & Laboratory Narrative & Table | 1 | 4 |
| 2. Sample Results and Raw Data (Organized By Sample) | 5 | 8 |
| a. ATL Sample Results Form | | |
| b. Target Compound Raw Data | | |
| -Internal Standard Area and Retention Time Summary (If Applicable) | | |
| -Surrogate Recovery Summary (If Applicable) | | |
| -Chromatogram(s) and Ion Profiles (If Applicable) | | |
| 3. QC Results and Raw Data | | |
| a. Method Blank (Results + Raw Data) | - | - |
| b. Surrogate Recovery Summary Form (If Applicable) | - | - |
| c. Internal Standard Summary Form (If Applicable) | - | - |
| d. Duplicate Results Summary Sheet | - | - |
| e. Matrix Spike/Matrix Spike Duplicate (Results + Raw Data) | - | - |
| f. Initial Calibration Data (Summary Sheet + Raw Data) | | |
| g. MDL Study (If Applicable) | - | - |
| h. Continuing Calibration Verification Data | | |
| i. Second Source LCS (Summary + Raw Data) | - | - |
| j. Extraction Logs | - | - |
| k. Instrument Run Logs/Software Verification | 9 | 13 |
| l. GC/MS Tune (Results + Raw Data) | - | - |
| 4. Shipping/Receiving Documents: | | |
| a. Login Receipt Summary Sheet | 14 | 15 |
| b. Chain-of-Custody Records | 16 | 17 |
| c. Sample Log-In Sheet | 18 | 19 |
| d. Misc. Shipping/Receiving Records (list individual records) | | |
| <u>Sample Receipt Discrepancy Report</u> | 20 | 22 |
| 5. Other Records (describe or list) | | |
| a. <u>Manual Spectral Defense</u> | - | - |
| b. <u>Manual Intergrations</u> | - | - |
| c. <u>Manual Calculations</u> | - | - |
| d. <u>Canister Dilution Factors</u> | - | - |
| e. <u>Laboratory Corrective Action Request</u> | - | - |
| f. <u>CAS Number Reference</u> | 23 | 24 |
| g. <u>Variance Table</u> | - | - |
| h. <u>Canister Certification</u> | - | - |
| i. <u>Data Review Check Sheet</u> | 25 | 25 |

Completed by:



(Signature)

Kara McKiernan/ Document Control

(Print Name & Title)

09/17/09

(Date)

WORK ORDER #: 0908455B

Work Order Summary

| | | | |
|------------------------|--|------------------|--|
| CLIENT: | Mr. Taeko Minegishi Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494 | BILL TO: | Accounts Payable Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494 |
| PHONE: | 800-825-5343 | P.O. # | 16512 |
| FAX: | 781-247-4305 | PROJECT # | 16512 |
| DATE RECEIVED: | 08/21/2009 | CONTACT: | Ausha Scott |
| DATE COMPLETED: | 09/16/2009 | | |

| <u>FRACTION #</u> | <u>NAME</u> | <u>TEST</u> |
|-------------------|----------------------|------------------|
| 21A | 100840 | ATL Applications |
| 21AA | 100840 Lab Duplicate | ATL Applications |
| 22A | 100839 | ATL Applications |
| 23A | 100838 | ATL Applications |
| 24A | 100177 | ATL Applications |
| 25A | 100178 | ATL Applications |
| 26A | 100180 | ATL Applications |
| 27A | 100181 | ATL Applications |
| 28A | 100331 | ATL Applications |
| 29A | 100332 | ATL Applications |
| 30A | 100333 | ATL Applications |
| 31A | 100334 | ATL Applications |
| 32A | 100335 | ATL Applications |
| 32AA | 100335 Lab Duplicate | ATL Applications |
| 33A | 100336 | ATL Applications |
| 34A | 100680 | ATL Applications |
| 35A | 100681 | ATL Applications |

Continued on next page


WORK ORDER #: 0908455B

Work Order Summary

| | | | |
|------------------------|--|------------------|--|
| CLIENT: | Mr. Taeko Minegishi Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494 | BILL TO: | Accounts Payable Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494 |
| PHONE: | 800-825-5343 | P.O. # | 16512 |
| FAX: | 781-247-4305 | PROJECT # | 16512 |
| DATE RECEIVED: | 08/21/2009 | CONTACT: | Ausha Scott |
| DATE COMPLETED: | 09/16/2009 | | |

| <u>FRACTION #</u> | <u>NAME</u> | <u>TEST</u> |
|-------------------|--------------|------------------|
| 36A | 100682 | ATL Applications |
| 37A | 100683 | ATL Applications |
| 38A | 100684 | ATL Applications |
| 39A | 100685 | ATL Applications |
| 40A | 100613 | ATL Applications |
| 41A | Method Blank | ATL Applications |
| 41B | Method Blank | ATL Applications |
| 41C | Method Blank | ATL Applications |
| 42A | CCV | ATL Applications |

CERTIFIED BY:



Laboratory Director

DATE: 09/16/09

This report shall not be reproduced, except in full, without the written approval of Air Toxics Ltd.

180 BLUE RAVINE ROAD, SUITE B FOLSOM, CA - 95630
(916) 985-1000 . (800) 985-5955 . FAX (916) 985-1020

LABORATORY NARRATIVE
Ozone by Radiello 172
Environmental Health & Engineering, Inc.
Workorder# 0908455B

Twenty Radiello 172 (Ozone) samples were received on August 21, 2009. The procedure involves reaction of 4-pyridylaldehyde with 3-methyl-2-benzothiazolinone hydrazone to yield the corresponding azide. The absorbance is then measured at 430 nm using a spectrophotometer. Results are reported in uG and uG/m3.

Sampling rate of 24.6 mL/min was provided by the manufacturer.

Receiving Notes

A Temperature Blank was included with the shipment. Temperature was measured and was not within 4 ± 2 °C. Coolant in the form of ice was present. Analysis proceeded.

Sample collection date was not provided on the Chain of Custody for all samples. The client was contacted and a dates were provided.

The Chain of Custody (COC) information for sample 100539 did not match the entry on the sample tag with regard to sample identification. The information on the COC was used to process and report the sample.

Analytical Notes

Results were calculated based on 25 deg C without temperature correction. The actual exposure time was used to calculate sample concentrations and reporting limits.

An exposure time of 22,000 minutes was used for the QC samples and samples 100336 and 100685.

All media used for the sampling were supplied by the client. Blank subtraction was not performed on the sample results since the media used for Method Blanks may be from a different lot than the media used for the samples.

Definition of Data Qualifying Flags

Eight qualifiers may have been used on the data analysis sheets and indicate as follows:

- B - Compound present in laboratory blank greater than reporting limit.
- J - Estimated value.
- E - Exceeds instrument calibration range.
- S - Saturated peak.
- Q - Exceeds quality control limits.
- U - Compound analyzed for but not detected above the detection limit.
- M - Reported value may be biased due to apparent matrix interferences.
- N - The identification is based on presumptive evidence.

File extensions may have been used on the data analysis sheets and indicates as follows:

a-File was requantified

b-File was quantified by a second column and detector

r1-File was requantified for the purpose of reissue

Sample Results and Raw Data

AIR TOXICS LTD.
ATL Application # 62 for RAD 172 (Ozone)
 Spectrophotometer

| Field Sample ID. | Lab Sample ID. | Collection Date | Analysis Date | Dilution Factor | Reporting Limit (ug) | Reporting Limit (ug/m3) | Amount (ug) | Amount (ug/m3) |
|----------------------|-------------------|--------------------|------------------|--------------------|-------------------------|----------------------------|----------------|-------------------|
| 100840 | 09084558-21A | 8/18/2009 | 8/21/2009 | 1.00 | 0.64 | 1.6 | 9.5 | 24 |
| 100840 Lab Duplicate | 09084558-21AA | 8/18/2009 | 8/21/2009 | 1.00 | 0.64 | 1.6 | 9.5 | 24 |
| 100839 | 09084558-22A | 8/18/2009 | 8/21/2009 | 1.00 | 0.64 | 1.6 | ND | ND |
| 100838 | 09084558-23A | 8/18/2009 | 8/21/2009 | 1.00 | 0.64 | 1.6 | ND | ND |
| 100177 | 09084558-24A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.2 | ND | ND |
| 100178 | 09084558-25A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.2 | ND | ND |
| 100180 | 09084558-26A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.2 | ND | ND |
| 100181 | 09084558-27A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.2 | ND | ND |
| 100331 | 09084558-28A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.5 | ND | ND |
| 100332 | 09084558-29A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.5 | ND | ND |
| 100333 | 09084558-30A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.5 | ND | ND |
| 100334 | 09084558-31A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.5 | ND | ND |
| 100335 | 09084558-32A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.5 | 9.0 | 21 |
| 100335 Lab Duplicate | 09084558-32AA | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.5 | 9.1 | 21 |
| 100336 | 09084558-33A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.2 | ND | ND |
| 100680 | 09084558-34A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.4 | ND | ND |
| 100681 | 09084558-35A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.4 | ND | ND |
| 100682 | 09084558-36A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.4 | 11 | 23 |
| 100683 | 09084558-37A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.4 | ND | ND |
| 100684 | 09084558-38A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.4 | ND | ND |
| 100685 | 09084558-39A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.2 | ND | ND |
| 100613 | 09084558-40A | 8/19/2009 | 8/21/2009 | 1.00 | 0.64 | 1.3 | ND | ND |
| Method Blank | 09084558-41A | NA | 8/21/2009 | 1.00 | 0.64 | 1.2 | ND | ND |
| Method Blank | 09084558-41B | NA | 8/21/2009 | 1.00 | 0.64 | 1.2 | ND | ND |
| Method Blank | 09084558-41C | NA | 8/21/2009 | 1.00 | 0.64 | 1.2 | ND | ND |
| CCV | 09084558-42A | NA | 8/21/2009 | 1.00 | 0.64 | 1.2 | %Rec 108 | |

COMMENTS: 1. NA=Not Applicable
 2. ND=Not Detected
 3. Exposure time of 21842 minutes was assumed for the QC samples and samples 100336 and 100685.
 4. Background subtraction not performed.

Ozone Radiello Calculation Worksheet

Workorder #: 09084558

Sampling Rate (ml/min): 24.6

Sampling T (deg C): 25

Volume (ml): 5

Date of Analysis: 8/21/2009

Corrected Q: 24.6

ically 24.6 for Ozone

Typically 25

typically 5 for Ozone

aking into account Temp

(Abs-Y-int)DF

Slope

Conc (ug) x 1000000

Q x Duration

Low PointDF

RL (ug) x 1000000

Q x Duration

LabSampleID

Client

Date of Collection

Abs

Duration (min)

DF

Ozone Conc (ug)

Conc (ug/m3)

RL(ug)

RL (ug/m3)

Result (ug)

| | | | | | | | | | | |
|------|----------------------|-----------|-------|-------|------|--------------|--------|-------|-------|-------------|
| 21A | 100840 | 8/18/2009 | 1.265 | 15840 | 1.00 | 9.474461336 | 24.314 | 0.638 | 1.638 | 9.474461336 |
| 21AA | 100840 Lab Duplicate | 8/18/2009 | 1.268 | 15840 | 1.00 | 9.497610953 | 24.374 | 0.638 | 1.638 | 9.497610953 |
| 22A | 100839 | 8/18/2009 | 0.044 | 15840 | 1.00 | 0.05256736 | 0.135 | 0.638 | 1.638 | ND |
| 23A | 100838 | 8/18/2009 | 0.050 | 15840 | 1.00 | 0.098865594 | 0.254 | 0.638 | 1.638 | ND |
| 24A | 100177 | 8/19/2009 | 0.058 | 21842 | 1.00 | 0.160598905 | 0.299 | 0.638 | 1.188 | ND |
| 25A | 100178 | 8/19/2009 | 0.059 | 21842 | 1.00 | 0.168831544 | 0.313 | 0.638 | 1.188 | ND |
| 26A | 100180 | 8/19/2009 | 0.068 | 21842 | 1.00 | 0.237764294 | 0.443 | 0.638 | 1.188 | ND |
| 27A | 100181 | 8/19/2009 | 0.066 | 21842 | 1.00 | 0.222331216 | 0.414 | 0.638 | 1.188 | ND |
| 28A | 100331 | 8/19/2009 | 0.045 | 17302 | 1.00 | 0.060283899 | 0.142 | 0.638 | 1.500 | ND |
| 29A | 100332 | 8/19/2009 | 0.046 | 17302 | 1.00 | 0.06800438 | 0.160 | 0.638 | 1.500 | ND |
| 30A | 100333 | 8/19/2009 | 0.044 | 17302 | 1.00 | 0.05256736 | 0.124 | 0.638 | 1.500 | ND |
| 31A | 100334 | 8/19/2009 | 0.048 | 17302 | 1.00 | 0.08343516 | 0.196 | 0.638 | 1.500 | ND |
| 32A | 100335 | 8/19/2009 | 1.209 | 17302 | 1.00 | 9.042335159 | 21.245 | 0.638 | 1.500 | 9.042335159 |
| 32AA | 100335 Lab Duplicate | 8/19/2009 | 1.210 | 17302 | 1.00 | 9.050051698 | 21.263 | 0.638 | 1.500 | 9.050051698 |
| 33A | 100336 | 8/19/2009 | 0.060 | 21842 | 1.00 | 0.176031982 | 0.328 | 0.638 | 1.188 | ND |
| 34A | 100880 | 8/19/2009 | 0.063 | 18412 | 1.00 | 0.19146506 | 0.440 | 0.638 | 1.409 | ND |
| 35A | 100881 | 8/19/2009 | 0.062 | 18412 | 1.00 | 0.19146506 | 0.440 | 0.638 | 1.409 | ND |
| 36A | 100882 | 8/19/2009 | 1.412 | 18412 | 1.00 | 10.60879255 | 23.422 | 0.638 | 1.409 | 10.60879255 |
| 37A | 100683 | 8/19/2009 | 0.053 | 18412 | 1.00 | 0.12201621 | 0.289 | 0.638 | 1.409 | ND |
| 38A | 100684 | 8/19/2009 | 0.073 | 18412 | 1.00 | 0.276346988 | 0.610 | 0.638 | 1.188 | ND |
| 39A | 100685 | 8/19/2009 | 0.037 | 21842 | 1.00 | -0.001448412 | -0.003 | 0.638 | 1.188 | ND |
| 40A | 100613 | 8/19/2009 | 0.061 | 19850 | 1.00 | 0.183748521 | 0.376 | 0.638 | 1.307 | ND |
| 41A | Method Blank | NA | 0.027 | 21842 | 1.00 | -0.078613801 | -0.146 | 0.638 | 1.188 | ND |
| 41B | Method Blank | NA | 0.029 | 21842 | 1.00 | -0.063180723 | -0.118 | 0.638 | 1.188 | ND |
| 41C | Method Blank | NA | 0.031 | 21842 | 1.00 | -0.047747645 | -0.089 | 0.638 | 1.188 | ND |
| 42A | CCV | NA | 0.933 | 21842 | 1.00 | 6.912570427 | 12.865 | 0.638 | 1.188 | 6.912570427 |

QC Duration

21842

CCV Spike Amt

6.384

Date of Calibration
8/21/2009 **Linear Regression**

0.129591779
0.037187702
0.998198851

QC Results and Raw Data

Spectrophotometer Logbook

@Air Toxics Ltd.

Log Book #: 1564Work Order: 0908455 B/CDate: 8/21/09Method: Rad 172Analyst: A. ToyamaWavelength: 430 nm

Prep. Notes:

| Standard ID | Concentration | ABS |
|------------------|-----------------------|-------|
| 1858 - 15 - 2.85 | 2.85 $\mu\text{g/mL}$ | 0.057 |
| 5.7 | 5.7 | 0.101 |
| 11.4 | 11.4 | 0.190 |
| 22.8 | 22.8 | 0.377 |
| 57 | 57 | 0.909 |
| 114 | 114 | 1.670 |

$$r = 0.99819$$

$$m = 0.12959$$

$$b = 0.037187$$

| Fraction | Dilution | ABS | Sample ID | Sample Volume |
|----------|----------|-------|-----------|---------------|
| 21A | 1.00 | 1.265 | 100840 | 5.0 mL |
| 22A | | 0.044 | 839 | |
| 23A | | 0.050 | 838 | |
| 24A | | 0.058 | 177 | |
| 25A | | 0.059 | 178 | |
| 26A | | 0.068 | 180 | |
| 27A | | 0.066 | 181 | |
| 28A | | 0.045 | 331 | |
| 29A | | 0.046 | 332 | |
| 30A | | 0.044 | 333 | |
| 31A | | 0.048 | 334 | |
| 32A | | 1.209 | 335 | |
| 33A | | 0.060 | 336 | |
| 34A | | 0.063 | 690 | |
| 35A | | 0.062 | 681 | |

Notes: _____

Spectrophotometer Logbook

@Air Toxics Ltd.

Log Book #: 1564Work Order: 0908455B/CDate: 8/21/09Method: Rad 172Analyst: A. ToyamaWavelength: 430 nmPrep. Notes: Cont. From page 12

| Standard ID | Concentration | ABS |
|--------------|-----------------------|-------|
| 1858-15-2.85 | 2.85 $\mu\text{g/mL}$ | 0.057 |
| 5.7 | 5.7 | 0.101 |
| 11.4 | 11.4 | 0.190 |
| 22.8 | 22.8 | 0.377 |
| 57 | 57 | 0.909 |
| 114 | 114 | 1.670 |

r = _____
 m = _____
 b = _____

| Fraction | Dilution | ABS | Sample ID | Sample Volume |
|----------|----------|---------|-----------|---------------|
| 36A | 1.00 | 1.412 | 100692 | 5.0 mL |
| 37A | | 0.053 | 683 | |
| 38A | | 0.073 | 684 | |
| 39A | | 0.037 | 685 | |
| 40A | | 0.061 | 613 | |
| 21A A | | 1.268 | 840 | |
| 32A A | | 1.210 | 335 | |
| BLK | | 0.027 | NA | |
| BLK | | 0.029 | | |
| BLK | | 0.031 | | |
| LCS/CCV | | 0.933 | | |
| | | 8/21/09 | | |
| | | AT | | |

Notes: Code 172 Lot 09146 Exp 01/10 used for BlanksLCS/CCV prepared at 57 $\mu\text{g/mL}$

Standard ID: 1858-14Project: Rad 172 MBTH SolutionAnalyst: A. ToyamaPreparation Date: 8/2/09Expiration Date: 8/2/09Solvent: DI H₂OSolvent Lot #: NAOzone

Procedure/Comments: Dissolve 2.5g of 3-Methyl-2-benzothiazolinone (Located
ERIA) hydrazine hydrochloride hydrate, 97% (1476-1106) into 500 mL DI
H₂O and add 2.5 mL of concentrated sulfuric acid.

8/20/09
AT

Standard ID: 1858-15

Project: Rad 172 Calibration Solution

Analyst: A. Toyama

Preparation Date: 8/21/09

Expiration Date: 8/21/09

Solvent: DI H₂O

Solvent Lot #: NA

Procedure/Comments: Dissolve 20 μ l of 4-Pyridine-carboxaldehyde, 97% (1476-1103, Located F22H) in 200 mL DI H₂O. From this solution prepare dilutions at 1:2, 1:5, 1:10, 1:20 and 1:40. = 114 μ g/mL

1:2) 250 μ l Pyridine solution with 250 μ l of DI H₂O. = 57 μ g/mL ✓

1:5) 100 μ l of Pyridine solution with 400 μ l of DI H₂O. = 22.8 μ g/mL

1:10) 100 μ l of Pyridine solution with 900 μ l of DI H₂O. = 11.4 μ g/mL

1:20) 250 μ l of Pyridine 1:10 solution with 250 μ l of DI H₂O. = 5.7 μ g/mL

1:40) 125 μ l of Pyridine 1:10 solution with 375 μ l of DI H₂O. = 2.85 μ g/mL

Then add 4.5 mL of MBTH solution to each level, stir and let stand for 1 hour (cover with parafilm) Then read the absorbance at 430 nm.

1 μ g of 4-pyridylaldehyde = 0.224 μ g of ozone

8/21/09

A.T.

Shipping/ Receiving Documents

**180 Blue Ravine Road, Suite B
Folsom, CA 95630**

**Phone (916) 985-1000 FAX (916) 985-1020
Hours 8:00 A.M. to 6:00 P.M. Pacific**

COMPANY: Environmental Health & Engineering, Inc.
ATTENTION: Mr. Taeko Minegishi
FAX #: 781-247-4305
FROM: Sample Receiving
Workorder #: 0908455B
of pages (Including Cover): 4

9/17/2009

Thank you for selecting Air Toxics Ltd. We have received your samples and have found discrepancies. In order to expedite analysis and reporting, please review the attached information for accuracy.

Corrections can be faxed to **Ausha Scott at 916-985-1020.**

ATL will proceed with the analysis as specified on the Chain of Custody and Sample Login page.

We have found a discrepancy between the Chain of Custody (COC) and the sample tag. The sample labeled 100539 on the COC is labeled as 100533 on the sample tag. ATL will report the sample identification on the COC unless otherwise notified.

Your prompt response is appreciated.

FROM: Environmental Health and Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494-2725

TO: AIR TOXICS

Please send invoices to ATTN: Accounts Payable
Please send reports to ATTN: Data Coordinator

In all correspondence regarding this matter, please refer to EH&E Project # 16512

The cost of this analysis will be covered by EH&E Purchase Order # 16512

For EH & E Data Coordinator - URGENT DATA ☒

| SAMPLE ID | SAMPLE TYPE | ANALYTICAL METHOD/NUMBER | OTHER:Time/Date/Vol. |
|-----------|-------------|--------------------------|----------------------|
| 100539 | AIR/PASSIVE | OZONE ANALYSIS | Ø |
| 100843 | | | 12 DAYS 11 DAYS |
| 100842 | | | |
| 100841 | | | |
| 100840 | | | |
| 100839 | | | |
| 100838 | | | |
| 100177 | | | 15D 4H 2MIN |
| 100178 | | | |
| 100180 | | | |
| 100181 | | | |
| 100331 | | | 12D 22MIN |
| 100332 | | | |
| 100333 | | | |
| 100334 | | | |
| 100335 | | | |

Special Instructions:

☒ Standard turn around time ☐ Rush by _____ date/time ☐ Other _____

☐ Fax results 781-247-4305

☐ RETURN SAMPLES

☒ Electronic transfer - datacoordinator@eh&e.com

☒ Additional report recipient mfragala@chem.com

Each signatory please return one copy of this form to the above address

Relinquished by: [Signature] of Environmental Health & Engineering, Inc. Date: 8/20/09

Received by: [Signature] 0830 of (company name) ATL Date: 8/21/09

Relinquished by: _____ of (company name) _____ Date: _____

Received by: _____ of (company name) _____ Date: _____

Relinquished by: _____ of (company name) _____ Date: _____

Received by: _____ of (company name) _____ Date: _____

Lab Data

Received by: _____ of Environmental Health & Engineering, Inc. Date: _____

Page 2 of 4

FROM: Environmental Health and Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494-2725

TO: AIR TOXICS

Please send invoices to ATTN: Accounts Payable
Please send reports to ATTN: Data Coordinator

In all correspondence regarding this matter, please refer to EH&E Project # 16512

The cost of this analysis will be covered by EH&E Purchase Order # 16512

For EH & E Data Coordinator - URGENT DATA ☒

| SAMPLE ID | SAMPLE TYPE | ANALYTICAL METHOD/NUMBER | OTHER:Time/Date/Vol. |
|------------|-------------|--------------------------|----------------------|
| 33A 100336 | AIR/PASSIVE | OZONE ANALYSIS | Ø |
| 34A 100680 | | | 12D 18H 52M |
| 35A 100681 | | | |
| 36A 100682 | | | |
| 37A 100683 | | | |
| 38A 100684 | | | |
| 39A 100685 | | | Ø |
| 40A 100613 | | | 13D 18H 50M |
| 41A 100614 | | | |
| 42A 100615 | | | |
| 43A 100616 | | | |
| 44A 100617 | | | |
| 45A 100618 | | | Ø |
| 46A 100619 | | | Ø |
| 47A 100165 | | | 13D 20H 49M |
| 48A 100166 | | | 1 |

Special Instructions:

☒ Standard turn around time

☐ Rush by _____ date/time


☐ Fax results 781-247-4305

☐ RETURN SAMPLES

☒ Additional report recipient

☒ Electronic transfer - datacoordinator@eh&e.com

mfragala@eh&e.com



Each signatory please return one copy of this form to the above address

Relinquished by: [Signature] of Environmental Health & Engineering, Inc. Date: 8/21/09

Received by: AD 0850 of (company name) ATC Date: 8/21/09

Relinquished by: _____ of (company name) _____ Date: _____

Received by: _____ of (company name) _____ Date: _____

Relinquished by: _____ of (company name) _____ Date: _____

Received by: _____ of (company name) _____ Date: _____

Lab Data

Received by: _____ of Environmental Health & Engineering, Inc. Date: _____

Page 3 of 4

SAMPLE RECEIPT SUMMARY

WORKORDER 0908455B

Client

Mr. Taeko Minegishi
Environmental Health &
Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494

Phone

800-825-5343

Fax

781-247-4305

Date Promised: 09/01/09 11:59 pm

Date Completed: 9/16/09

Date Received: 8/21/09

PO#: 16512

Project#: 16512

Total \$: \$ 1,100.00

Sales Rep: TL

Logged By: MG

| <u>Fraction</u> | <u>Sample #</u> | <u>Analysis</u> | <u>Collected</u> | <u>Amount\$</u> |
|-----------------|----------------------|------------------|------------------|-----------------|
| 21A | 100840 | ATL Applications | 8/18/2009 | \$50.00 |
| 21AA | 100840 Lab Duplicate | ATL Applications | 8/18/2009 | \$0.00 |
| 22A | 100839 | ATL Applications | 8/18/2009 | \$50.00 |
| 23A | 100838 | ATL Applications | 8/18/2009 | \$50.00 |
| 24A | 100177 | ATL Applications | 8/19/2009 | \$50.00 |
| 25A | 100178 | ATL Applications | 8/19/2009 | \$50.00 |
| 26A | 100180 | ATL Applications | 8/19/2009 | \$50.00 |
| 27A | 100181 | ATL Applications | 8/19/2009 | \$50.00 |
| 28A | 100331 | ATL Applications | 8/19/2009 | \$50.00 |
| 29A | 100332 | ATL Applications | 8/19/2009 | \$50.00 |
| 30A | 100333 | ATL Applications | 8/19/2009 | \$50.00 |
| 31A | 100334 | ATL Applications | 8/19/2009 | \$50.00 |
| 32A | 100335 | ATL Applications | 8/19/2009 | \$50.00 |
| 32AA | 100335 Lab Duplicate | ATL Applications | 8/19/2009 | \$0.00 |
| 33A | 100336 | ATL Applications | 8/19/2009 | \$50.00 |
| 34A | 100680 | ATL Applications | 8/19/2009 | \$50.00 |
| 35A | 100681 | ATL Applications | 8/19/2009 | \$50.00 |
| 36A | 100682 | ATL Applications | 8/19/2009 | \$50.00 |
| 37A | 100683 | ATL Applications | 8/19/2009 | \$50.00 |
| 38A | 100684 | ATL Applications | 8/19/2009 | \$50.00 |

Note: Samples received after 3 P.M. PST are considered to be received on the following work day.
Atlas Project Name/Profile#: CPSC Indoor Air Monitoring/13297

BILL TO: Accounts Payable
Environmental Health & Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494

Analysis Code: Other GC

TERMS:

Reporting Method: ATL Application #62 Ozone-Radiello 172

180 BLUE RAVINE ROAD, SUITE B FOLSOM, CA - 95630
(916) 985-1000 . (800) 985-5955 . FAX (916) 985-1020

SAMPLE RECEIPT SUMMARY Continued

| | | |
|------------------------|--------------|---|
| Client | Phone | Date Promised: 09/01/09 11:59 pm |
| Mr. Taeko Minegishi | | Date Completed: 9/16/09 |
| Environmental Health & | 800-825-5343 | Date Received: 8/21/09 |
| Engineering, Inc. | Fax | PO#: 16512 |
| 117 Fourth Avenue | | Project#: 16512 |
| Needham, MA 02494 | 781-247-4305 | |
| Sales Rep: TL | | Total \$: \$ 1,100.00 |
| | | Logged By: MG |

| <u>Fraction</u> | <u>Sample #</u> | <u>Analysis</u> | <u>Collected</u> | <u>Amount\$</u> |
|--|-----------------|------------------|------------------|-----------------|
| 39A | 100685 | ATL Applications | 8/19/2009 | \$50.00 |
| 40A | 100613 | ATL Applications | 8/19/2009 | \$50.00 |
| 41A | Method Blank | ATL Applications | NA | \$0.00 |
| 41B | Method Blank | ATL Applications | NA | \$0.00 |
| 41C | Method Blank | ATL Applications | NA | \$0.00 |
| 42A | CCV | ATL Applications | NA | \$0.00 |
| Misc. Charges eCVP (20) @ \$5.00 each. | | | | \$100.00 |

Note: Samples received after 3 P.M. PST are considered to be received on the following work day.
Atlas Project Name/Profile#: CPSC Indoor Air Monitoring/13297

BILL TO: Accounts Payable
Environmental Health & Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494

Analysis Code: Other GC

TERMS:

Reporting Method: ATL Application #62 Ozone-Radiello 172

180 BLUE RAVINE ROAD, SUITE B FOLSOM, CA - 95630
(916) 985-1000 . (800) 985-5955 . FAX (916) 985-1020

Sample Discrepancy Report

Identification

Initiated By: MG Project ID: 13297 PM: BL Date: 8/21/2009 Discrepancy Type: ☒ 1. ☒ 2. ☐ 3.

Workorder(s) affected: 0908455 Sample(s) affected: ALL, 17A

1. Sample Receipt Discrepancies

Narration Not Required:

- 1.1. ☐ Sample container (cartridge/tube/VOA vial) was received broken, however sample was intact.
- 1.2. ☐ No brass cap on canister.
- 1.3. ☐ Date of Collection noted on first sample, but no arrow down to indicate all samples.

Notify Lab for further determination:

- 1.4. ☐ Tedlar bag received with minimal volume.

Initials: _____ Date: _____

Narration Required in Lab Narrative and Sample Confirmation:

- 1.5. ☐ COC was not filled out in ink.
- 1.6. ☐ COC improperly relinquished / received.
- 1.7. ☒ Sample tags / can numbers do not match the COC.
- 1.8. ☐ Sample date ☐ error / ☐ missing on COC but noted on sample tag (check one).
- 1.9. ☐ Custody Seal on the outside of the container was ☐ broken / ☐ improperly placed (check one).
- 1.10. ☐ ID-none on the sample Tag/Blank
- 1.11. ☐ Other (describe below).

Describe the Discrepancy: 1.7: Sample 17A: The Sample ID tag reads "100533"

2. Sample Receipt/Screening Discrepancies requiring PM notification

Document on Cover Page of Sample Receipt Confirmation and in Receiving Notes of Lab Narrative

If Section II. is filled out PM must be notified within 24 hrs of initiation

- 2.1. ☐ COC was not received with samples.
- 2.2. ☐ Analysis method(s) is ☐ not specified / ☐ incorrectly specified (check one) on the COC.
- 2.3. ☐ Incorrect sampling media / container for analysis requested.
- 2.4. ☐ Number of samples on the COC does not match the number of samples that were received.
- 2.5. ☐ Samples were received expired.
- 2.6. ☒ Sampling date (time for sulfur) is not documented for ☐ some / ☒ any samples (check one).
- 2.7. ☐ Sample received with amount of H₂O in the Tedlar Bag.
- 2.8. ☐ Sample cannot be analyzed. Container was ☐ received broken / ☐ leaking / ☐ flat / ☐ defective.
- 2.9. ☐ Tedlar bag / canister received emitting a strong odor; Sample ☐ can / ☐ cannot (check one) be analyzed.
- 2.10. ☐ Tedlar Bag for Sulfur analysis has metal fitting.
- 2.11. ☐ Environmental Supply Company valves
- 2.12. ☐ Sorbent samples-sampling volume was not provided
- 2.13. ☐ Flow controller used – canister samples received at ambient or under pressure.
- 2.14. ☐ Canister was at ambient pressure at time of pressurization and (check all that apply):
- ☐ Canister failed leak check on two manifolds,
 - ☐ Canister valve was open,
 - ☐ Brass nut was loose/not present.
 - ☐ Sample can be analyzed
 - ☐ Cannot be analyzed
- 2.15. ☐ Canister sample received with a vacuum difference >5.0"Hg between the receipt vac. And the final vac. reported on the COC, indicating loss of vacuum.
- 2.16. ☐ Canister sample received at >15"Hg (not identified as a Trip/Field Blank).
- 2.17. ☐ Canister Trip Blank received at low vacuum (< 25"Hg).
- 2.18. ☒ Sorbent Sample received outside method required temperature of 2°C to 6°C; ☒ ice / ☒ blue ice (check one) was present. A temp. Blank ☐ was / ☒ was not present (check one).
- 2.19. ☐ Other (describe below)

Initials: _____

Date: _____

Notify Receiving: ☐

Notify PM: ☐

Describe the Discrepancy: 2.6: Date of Collection not noted on the COC.

2.18: Samples arrived at 8C

3. Lab Discrepancies requiring Team Leader/PM notification

Document in Analytical Notes of Lab Narrative

If Section III. is filled out PM must be notified within 24 hrs of initiation

- | | |
|--|--|
| 3.1. <input type="checkbox"/> Tedlar Bag found to be leaking at the time of analysis; sample <input type="checkbox"/> can / <input type="checkbox"/> cannot (check one) be analyzed. | 3.6. <input type="checkbox"/> Sample loss due to instrument malfunction / broken glassware. |
| 3.2. <input type="checkbox"/> Tedlar Bag found to be flat/low volume; sample cannot be analyzed. | 3.7. <input type="checkbox"/> Low/high surrogate recoveries noted in QC/sample(s) for extractable samples. |
| 3.3. <input type="checkbox"/> Sulfur samples received with insufficient time to analyze prior to expiration. | 3.8. <input type="checkbox"/> Reporting Limit was raised. |
| 3.4. <input type="checkbox"/> Canister found to be leaking at the time of analysis. | 3.9. <input type="checkbox"/> Post weight > Pre weight in field/lab Blank for PM10/TSP samples. |
| 3.5. <input type="checkbox"/> VOST tube saturated; bag dilution necessary. | 3.10. <input type="checkbox"/> Other (describe below). |

Initials: _____ Date: _____ Notify Receiving: ☐ Notify PM: ☐

Team Lead Initials: _____ Date: _____

Describe the Discrepancy: _____

How Does this Affect Client: _____

Project Manager Use Only

Project Manager Notification Complete

☒ Section 2 Complete

☐ Section 3

Action:

- ☐ It is not necessary to notify the client. Narrate the discrepancy in Receiving Notes/Analytical Notes of Lab Narrative.

PM Initials: _____ Date: _____

- ☒ Client notification required. See attached client contact / email, or comments below:

Client Notification:

PM Initials: BL Person notified: David Shore

Date: 8/21/2009

- ☐ Waiting for Client Reply

Comments: **Proceed and narrate temperature discrepancy. See table for time of collection.**

☐ Notify Lab Name: _____ Date: _____ Notify Receiving: ☒

- ☐ Additional notifications attached.

Additional Comments:

Other Records

Method : ATL Application #62 Ozone-Radiello 172

| CAS Number | Compound | Rpt. Limit (ug) |
|------------|----------|-----------------|
| 10028-15-6 | Ozone | 1.0 |

DATA REVIEW CHECKLIST

Work Order #:

0908455B

| | | | | | | |
|--------------------------|--------------------------|-------------------------------------|-------------------------------------|-------------------------------------|--------------------------|--|
| A ₁ | A ₂ | R | T | M | Q | |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Analysis/Reporting vs. Project Profile/SOP requirements checked (i.e. 100% Dups, J-Flag to MDL, etc) |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | The final report has the correct reporting list, special units, and header info. |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Lab Narrative is correct (proper method & description/Receiving & Analytical notes correct) |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Sample Discrepancy Report (SDR) is completed |
| <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Corrective Action issued - # _____ |
| <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Unusual circumstances have been documented in the notes section below |

LUMEN validation report present and initialed

CIRCLE (YES (NO))

| | | | | | | |
|--------------------------|--------------------------|-------------------------------------|--------------------------|-------------------------------------|--------------------------|---|
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Lab Blank, CCV, LCS and DUP met QC criteria |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Hold time is met for all samples |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Appropriate data qualifier flags are applied |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Manual integrations for samples and QC are properly documented |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Samples analyzed within the project or method specific clock |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Retention times have been verified |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Appropriate ICAL(s) included |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | At least one result per sample is verified against the target quant sheets/raw data |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Dilution factor correctly calculated (sample load volume, syringe and bag dilutions, can pressurization(s)) |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Correct amount of sample analyzed (i.e. sample not over-diluted) |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Spectra verified - documentation of spectral defense included (Section 5A of eCVP pkg) |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | TICs resemble reference spectra |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | TICs between duplicate samples are consistent |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Checked samples for trends (i.e. Influent vs. Effluent, Field Dups, Field/Trip Blank, etc.) |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Data for multiple analyses of sample(s) has been evaluated for comparability of results |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Special units for all samples in the final report are correctly calculated |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Manually entered results checked (i.e. TPH/NMOC) |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Chain of Custody verified for any special comments (i.e. different compounds/RLs, action levels) |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Chain of Custody scanned correctly |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Verify sample id's vs. chain of custody |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Date MDL(s) performed per instrument(s) |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Samples pressurized w/ appropriate gas (N ₂ or He) <input checked="" type="checkbox"/> Other (i.e. Tedlar bag, cartridge, sorbent) |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Final pressure consistent with canister size (6L vs. 1L) |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Verify receipt pressures |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Verify canister ID #'s |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Final invoice amount correct (adjusted for TAT, Penalties, Re-issue Charges etc.) |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | MDL date(s) present for all instruments utilized |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Client LUMEN report reviewed for accuracy and completeness |

Notes: (to include: noting samples with QA/QC problems, Blanks with positive hits, narratives, etc.)

A/R:

Dup. 21A, 32A

M/Q:

| | | | |
|--------------------------------|-------------------------|--------------------------|------------------|
| A ₁ /A ₂ | R/T | M | Q |
| (Analytical Review/Date) | (Reporting Review/Date) | (Management Review/Date) | (QA Review/Date) |
| A ₁ : | R: 4/9/16/09 | M: 4/16/09 | Q: |
| A ₂ : | T: | | |

Note (1): Please check all the appropriate boxes. Indicate "NA" for any statement that does not apply.

Rev. 02/20/09

Note (2): Management reviewer and reporting reviewer must be separate individuals.